

| | Final body weight (g) | Minimum floor area (cm ²) | Floor area per animal (groups) | Minimum height (cm) |
|---------------------|-----------------------|---------------------------------------|--------------------------------|---------------------|
| EUROPE ¹ | < 600 | 4500 | 1500 | 50 |
| | > 600 | 4500 | 3000 | 50 |
| | adult males | 6000 | 6000 | 50 |
| | jill and litter | 5400 | 5400 | 50 |
| UK ² | < 800 | 2250 | 1500 | 50 |
| | > 800 | 4500 | 3000 | 50 |

1 Requirements for Establishments and for the Care of Accommodation of Animals, Annex III of the Directive 2010/63/EU of the European Parliament and of the Council on the protection of animals used for scientific purposes, 22 September 2010

2 Home Office - Animals (Scientific Procedures) Act 1986 Amendment Regulations 2012

Cage specifications and dimensions

"P" TYPE SOLID CAGE BASE

- Size (L x W x H): 870 x 740 x 130 mm
34.25 x 29.13 x 5.118 in
- Floor area: 5456 cm² / 5.87 ft²
- Internal Height: 600 mm / 23.62 in

Rack specifications and dimensions

- N. of Cages: 2 (1W x 2H)
- Overall Rack Dimensions (L x W x H):
940 x 845 (+160) x 1700 mm
37.01 x 33.27 (+6.299) x 66.93 in



ENVIRONMENTAL ENRICHMENT: RAISED SHELF



ENVIRONMENTAL ENRICHMENT: INTERLINKING TUBE



Animal Welfare

- The rack interlinking increases floor space and allows communication between cages for social housing
- The interlinking tube connects two cages to enhance social behaviour
- The raised shelf and nesting box, together with the internal height, create an ideal three-dimensional enriched environment
- The newly designed stainless-steel rack with noise-free Polycarbonate panels helps to avoid animal stress
- Transparent polycarbonate gate avoids air drafts on the animal

Ergonomics, Functionality and Cost-effectiveness

- A push-and-close system allows to easily close the gate even while holding animals; the gate features also a safe locking system to avoid animals escaping from the cage
- A fast-dry roof design results in quicker and easier washing procedures
- Rounded corners and lighter rack ensure easy and safe transportation.

Visibility: easy daily check

Transparent removable side panels allow optimal visibility of the entire cage volume: quick and easy daily check

