

EASY-BOX  
CO2 CHAMBER  
RACK SYSTEM



POST-  
TREATMENT




AUTOMATED SYSTEM FOR EUTHANIZING ANIMALS  
AFTER EXPERIMENTATION



# POST-TREATMENT

System for euthanising animals after experimentation in full compliance with legislation and ethical standards.



## OPERATION DETAILS

The processes of all our systems are **automated, safe** and **strictly reproducible**. They consist of a filling, a maintenance and a gas extraction phase.



Compliant with European Directive 2010/63/EU

## TACTILE CONTROL PANEL

The control panel can be used to select the most appropriate cycle for the **species to be treated**.

It allows the process to be monitored **in real time**. **Safety alarms** interrupt the cycle if an anomaly is detected.



## CONTROL & VERIFICATION

For all our systems, the concentrations of CO<sub>2</sub> are **measured** and **calibrated** in the factory.

A **specific verification** report is provided at the delivery time and after every on-site maintenance operation.



### RACK

This modular CO2 chamber system offers an adaptable **treatment volume**.

Each chamber is **independent** and has its own cycle controller. **Large indicators** allow you to see the availability of each chamber at a glance. These units can be placed on a bench or mounted onto a **mobile three-place stand**.

### MC2003 CO2 CHAMBERS

Designed to be fully **hermetic**, these chambers come in several sizes, **respectively 38, 120 and 180 litres**. These independent chambers are designed to accept different **cages** in standard Type II or Type III formats. They are also available in a **mobile version** on a specially equipped trolley.



### EASY-BOX

This new system helps simplify euthanasia procedures by **limiting the stress experienced by the animal**. The process takes place **without removing animals** from their cage. The CO2 enrichment system operates **very silently**.

The system is very simple to use: simply place the cage on the table, remove the grill, position the **Easy-Box** over the cage and start the cycle.

# SIMPLE & LONG LASTING SOLUTIONS

A collaborative approach to help with  
your projects...as simple as that!



Z.I Rue de l'Aube - 51 310 ESTERNAY - FRANCE  
Tél. +33 (0)3 26 42 50 15 - fax +33 (0)3 26 42 50 16  
email : [contact@minerve.eu](mailto:contact@minerve.eu)

[www.minerve.eu](http://www.minerve.eu)

Equipement Vétérinaire MINERVE SA - RCS REIMS 324755545

